


Medium-Temperature and High-Pressure Glass-Based pH Probes

Our patent-pending glass-based medium-temperature and high-pressure pH electrodes are divided into three types (A, B, and C).

Type A or B high pressure pH electrode has a compression fitting with NPT* thread and an insertion tube. The fitting of **Type A** has a metal seal that will lock to the insertion tube upon first use (insertion depth adjustable one-time only). The fitting of **Type B** has a soft seal that can be repositioned on the insertion tube after each use (insertion depth fully adjustable).

Type C pH electrode has a thread fitting. If the thread size is 1/4" or larger, a thin-walled protection tube may be provided. Type C design is ultra compact and may be used in autoclaves that have smaller access ports (1/8" or 1/4" NPT). The insertion depth for Type C high pressure pH electrode is fixed at factory (non adjustable).


How To Order Medium-Temperature High-Pressure Glass-Based pH Electrodes Example



High-Temperature and High-Pressure Platinum ORP/Conductivity Probes

High-temperature and high-pressure platinum probes for redox potential (ORP) and conductivity measurements are divided into two types (A and B)

Type A or B platinum probe has a compression fitting with NPT* thread and an insertion tube. The fitting of **Type A** has a metal seal that will lock to the insertion tube upon first use (insertion depth one-time adjustable only). The fitting of **Type B** has a soft seal that can be repositioned on the insertion tube after each use (insertion depth fully adjustable).




*NPT - American National Pipe Thread.

High-Temperature & High-Pressure pH, ORP and Reference Probes

Reference Electrodes, pH Electrodes, and Platinum Electrodes for Redox Potential (ORP) and Conductivity Measurements




Corr Instruments, LLC
303 Clarence Tinker Dr.
San Antonio, TX 78226, USA



Phone: 210 748-4073
Fax: 210 696-1126
Email: info@corrinstruments.com
Web: www.Corrinstruments.com

Corr Instruments





Corr Instruments offers a full range of exceptionally high performance, high-temperature and high-pressure pH electrodes, reference electrodes, platinum electrodes for redox potential (ORP) and conductivity measurements at temperatures up to 650°F (343°C) and pressures up to 5100 psi (35 MPa). These probes are ideal tools for corrosion monitoring and thermodynamic/electrochemical studies in harsh environments such as steam generators, nuclear power reactors, geothermal systems, deep oil/gas wells, and super critical water systems.

Many of these probes are built on our patent-pending technology and our innovative sealing compound, Queon™. Queon is the only known material in the world for sealing electrical conductors in metal sealing glands or compression feedthroughs that operate at temperatures up to 650°F (343°C).



Probe Selection Guide

Image	Probe Type	Temperature Rating	Pressure Rating
	Ag/AgCl Reference Probe	0 to 305 °C (32 to 581 °F)	up to 5100 psi (350 atm)
	ZrO ₂ -Based pH Probe	130 to 305 °C (366 to 650 °F)	2000 Psi (136 atm)
	Glass-Based pH Probe	0 to 120 °C (32 to 248 °F)	2000 psi (136 atm)
	Platinum Redox (ORP) Probe	0 to 305 °C (32 to 581 °F)	up to 5100 psi (350 atm)
	Platinum Conductivity Probes	0 to 305 °C (32 to 581 °F)	up to 5100 psi (350 atm)

Product Range and Selection

High-Temperature and High-Pressure Pressure-Balanced Ag/AgCl Reference Probes


Our patent-pending high-temperature/high-pressure Pressure-Balanced Ag/AgCl Reference Electrodes are divided into three types (A, B, and C) and two designs (External Reference Electrode and Internal Reference Electrode)

Type A or B has a compression fitting with NPT* thread and an insertion tube. The fitting of **Type A** has a metal seal that will lock to the insertion tube upon first use (insertion depth adjustable one-time only). The fitting in **Type B** has a soft seal that can be repositioned on the insertion tube after each use (insertion depth fully adjustable).

Type C has a threaded fitting. If the thread size is 1/4" or larger, a thin-walled protection tube may be provided. Type C is ultra compact and may be used in autoclaves that have smaller access ports (1/8" NPT). The insertion depth for Type C is fixed at factory (non adjustable). Type C is available for internal design only.

External design: Reference material (Ag/AgCl) is away from heated location and maintained near ambient temperature so that the probe is more stable when used at T>200 °C.

Internal design: Reference material (Ag/AgCl) is in heated location so that the probe has a better defined thermodynamic potential. Not recommended for long-term use at T>200 °C.




High-Temperature and High-Pressure Zirconia-Based pH Probes

ZrO₂-based high-temperature and high-pressure pH electrodes are divided into three types (A, B, and C)

Type A or B pH electrode has a compression fitting with NPT* thread and an insertion tube. The fitting of **Type A** has a metal seal that will lock to the insertion tube upon first use (insertion depth adjustable one-time only). The fitting in **Type B** has a soft seal that can be repositioned on the insertion tube after each use (insertion depth fully adjustable).

Type C pH electrode has a thread fitting. If the size is 1/4" or larger, a thin-walled protection tube may be provided. Type C design is ultra compact and may be used in autoclaves that have smaller access ports (1/8" or 1/4" NPT). The insertion depth for Type C is fixed at factory (non adjustable).



*NPT -American National Pipe Thread.